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General Information about the method

After alkaline decomposition tin is measured using atomic emission spectrometry. The usual method to decompose H₂O₂ in a platinum dish cannot be used here because tin sticks to the platinum thus reducing the amount found in the analysis. For the tin determination an alkaline decomposition is therefore carried out in a glass beaker. With the help of an atomizer the sample solution is atomized and is then transported in the form of an aerosol using a carrier gas (argon = atomizing gas) into the inductively coupled plasma (ICP). In the plasma the sample is vapourized, the molecules split into atoms and the atoms excited to emission. Due to the high plasma temperature (approx. 8000°C) ionic lines as well as atomic lines occur and are also used in the analysis. The emission lines are separated spectrally in the spectrometer. The intensities of the emitted element lines are measured with a photomultiplier.

Equipment

- atomic emission spectrometer
- microlitre-pipettes (Eppendorf/Brand)
- heating plate

Reagents

- hydrogen peroxide
- dest. Water - prepared in quartz distillation apparatus
- ICP standard solutions - ready to use standard solutions available from different suppliers
- hydrochloric acid 30 % SUPRAPUR®
- sodium hydroxide solution c(NaOH) ~ 5 mol/l
- argon 4.8 (99.998 %)

Analytical Methods for Hydrogen Peroxide

Determination of tin using atomic emission spectrometry (ICP-AES)

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Special safety instructions

All reagents and chemicals must be handled according to the health and safety regulations. Refer to the safety data sheets.

Special procedure instructions

Danger of decomposition by contact with incompatible materials, contaminants, metals, reducing agents, alkalis.

Procedure

Sample preparation:

Place 2.5 ml sodium hydroxide solution c(NaOH) ~ 5 mol/l into a scratch-free 1 l beaker and carefully add 25 ml H₂O₂ in small portions. In the beginning warm until the decomposition reaction starts. Should the decomposition become too violent add some distilled water (wash bottle) to cool it down. After all the H₂O₂ has been decomposed acidify with 2.5 ml hydrochloric acid 0% and rinse into a 50 ml volumetric flask and fill to the mark.

Calibration:

For the calibration a blank and a standard solution containing the same reagents (sodium hydroxide solution and hydrochloric acid) in the same quantities as in the samples are prepared. The standard solution also contains 5mg/l Sn for the calibration. The calibration is checked using a check solution of 1 mg/l Sn and the same matrix as described above. Before calibrating always carry out an "Auto Search". As a result of temperature differences and mechanical faults, the wavelengths found with the spectrometer do not correspond exactly to the literature values. "Auto Search" stores the deviations (Offsets), so that the literature wavelengths can be used.

Analytical Methods for Hydrogen Peroxide

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Measurement:

Using the coefficients determined during the calibration the intensities measured can be converted to concentrations.

Should the concentration be higher than that of the calibration standards, repeat the measurement using a suitable calibration

Calculation

$c(\text{real}) \text{ (mg/l)} = c(\text{measured}) * \text{dilution factor}$

$$c(\text{real}) \text{ (ppm)} = \frac{c(\text{measured}) * \text{dilution factor}}{\text{density}}$$

accordance with local regulations

Literature

- Manufacturer equipment description
- Product information "Hydrogen Peroxide"
- Joachim Nölte, ICP Emissionsspektrometrie für Praktiker, Wiley-VCH Verlag GmbH

Remarks

The method is based on the internal analytical method WM21.

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.